

EPSA[®] *safety first*

**CC TRACHEOSTOMY
CANNULA WITH CUFF**

CC TRACHEOSTOMY CANNULA WITH CUFF

CODE: 9755

1 PVC cannula, conical shape, curved, longitudinal line opaque to X-rays, polished distal orifice, atraumatic. Pivoting fixation flange. Low pressure cuff, pilot balloon and retention valve.

Pilot balloon and retention valve

2 inner PVC cannulas, interchangeable

1 insertion aid

1 cough cap

1 standard connector.

1 velcro neck band



Available Sizes

5,0 / 5,5 / 6,0 / 6,5 / 7,0 / 7,5 / 8,0 /
8,5 / 9,0 / 9,5 / 10,0 mm

STERILE, SINGLE PACK, SINGLE USE.

Recommendation for use

Recommendation of a professional on the use of the technology.

The CC Tracheostomy Cannula with cuff is indicated in patients who require mechanical ventilation due to functional, post-traumatic or surgical alterations that affect the patency of this sector of the airway.

Its purpose is:

- enable mechanical ventilation;
- allow endoscopic maneuvers (diagnostic and / or therapeutic in the airway);
- facilitate the humidification of secretions and endotracheal aspiration.

The cuff seals the tracheal lumen, allowing adequate control of gas movement during mechanical ventilation and preventing the passage of fluids and secretions from the upper airway into the bronchial tree.



Recommendation for use

The product is offered with:

- (2) inner cannulas, so that the patient always has one in use and another clean for their substitution,
 - (1) standard connector, which allows the connection of the tracheostomy tube to:
 - manual or automatic mechanical ventilation devices (respirators),
 - nebulization devices.
 - (1) cough cap, used to prevent tracheo-bronchial secretions from being ejected during expiration or coughing spells.
 - (1) velcro neck band (velour with velcro), made of material whose texture ensures a correct fixation of the cannula without risk of breaking the plate.
- It is supplied with an introduction mandrel to facilitate its application in recently performed tracheostomies.

