



EPSA[®]


ELECTROPLAST S.A.

CATALOGUE



EPSA is a leading Uruguayan company in the manufacturing of medical devices. The company was founded in 1960 and is presently the most important local supplier in its field. Over the years, EPSA has gained prestige and acknowledgement while satisfying the clients' needs and responding to agrowing demand. Using the most advanced technology and delivering superior quality products allowed EPSA to gain access to regional and international markets. The company is based in Montevideo-Uruguay and has a 3.000 m2 production plant, including a modern Ethylene Oxide sterilization chamber, administrative offices and a warehouse. The company currently employs about 130 people working under strict hygiene standards. EPSA has GMP, EN ISO13485 and EC Certificates.

CERTIFICATES



TRACHEOSTOMY CANNULA WITHOUT CUFF

PRODUCT CODE: 9751



- ◆ 1 PVC cannula, curved, flexible, transparent, longitudinal line opaque to X-rays, polished tip, atraumatic.
Flexible fixation flange with standard connector.
- ◆ 1 insertion aid with olive tip.
- ◆ 1 velcro neck band.

STERILE, SINGLE PACK, SINGLE USE.

AVAILABLE SIZES: 5,0 / 5,5 / 6,0 / 6,5 / 7,0 / 7,5 / 8,0 / 8,5 / 9,0 / 9,5 / 10,0 mm.



PEDIATRIC AND NEONATAL TRACHEOSTOMY CANNULA WITHOUT CUFF

PRODUCT CODE: 9764 (PEDIATRIC) - 9758 (NEONATAL)

PEDIATRIC

- ◆ 1 PVC cannula, curved, flexible, transparent, longitudinal line opaque to X-rays, polished tip, atraumatic.
Flexible fixation flange with standard connector.
- ◆ 1 cloth tape for neck fixation.

STERILE, SINGLE PACK, SINGLE USE.

AVAILABLE SIZES: 3,0 / 3,5 / 4,0 / 4,5 mm.

NEONATAL

- ◆ 1 PVC cannula, curved, flexible, transparent, longitudinal line opaque to X-rays, polished tip, atraumatic.
Flexible fixation flange with standard connector.
- ◆ 1 cloth tape for neck fixation.

STERILE, SINGLE PACK, SINGLE USE.

AVAILABLE SIZES: 2,5 / 3,0 / 3,5 / 4,0 / 4,5 mm.



TRACHEOSTOMY CANNULA WITH CUFF

PRODUCT CODE: 9722



- ◆ 1 PVC cannula, curved, flexible, transparent, longitudinal line opaque to X-rays, polished tip, atraumatic.
Flexible fixation flange with standard connector.
Low-pressure cuff, pilot balloon and retention valve.
- ◆ 1 insertion aid with olive tip.
- ◆ 1 velcro neck band.

STERILE, SINGLE PACK, SINGLE USE.

AVAILABLE SIZES: 5,0 / 5,5 / 6,0 / 6,5 / 7,0 / 7,5 / 8,0 / 8,5 / 9,0 / 9,5 / 10,0 mm.



PEDIATRIC TRACHEOSTOMY CANNULA WITH CUFF.

PRODUCT CODE: 9766



- ◆ 1 PVC cannula, curved, flexible, transparent, longitudinal line opaque to X-rays, polished tip, atraumatic. Flexible fixation flange with standard connector. Low-pressure cuff, pilot balloon and retention valve.
- ◆ 1 cloth tape for neck fixation.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 3,0 / 3,5 / 4,0 / 4,5 mm.



CC TRACHEOSTOMY CANNULA WITHOUT CUFF

PRODUCT CODE: 9731

- ◆ 1 PVC cannula, conical shape, curved, polished distal orifice, atraumatic. Pivoting fixation flange.
- ◆ 2 inner PVC cannulas, interchangeable.
- ◆ 1 phonation valve.
- ◆ 1 cough cap.
- ◆ 1 standard connector.
- ◆ 1 velcro neck band.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 4,0 / 5,0 / 6,0 / 7,0 / 8,0 / 9,0 / 10,0 / 11,0 / 12,0 mm.



CC TRACHEOSTOMY CANNULA WITH INTRODUCER WITHOUT CUFF

PRODUCT CODE: 9735




- ◆ 1 PVC cannula, conical shape, curved, polished distal orifice, atraumatic. Pivoting fixation flange.
- ◆ 2 inner PVC cannulas, interchangeable.
- ◆ 1 insertion aid.
- ◆ 1 phonation valve.
- ◆ 1 cough cap.
- ◆ 1 standard connector.
- ◆ 1 velcro neck band.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 4,0 / 5,0 / 6,0 / 7,0 / 8,0 / 9,0 / 10,0 / 12,0 mm.



CC FENESTRATED TRACHEOSTOMY CANNULA WITHOUT CUFF

PRODUCT CODE: 9251



- ◆ 1 PVC cannula, conical shape, curved, polished distal orifice, atraumatic. Phonation opening. Pivoting fixation flange.
- ◆ 2 inner PVC cannulas, interchangeable, with phonation opening.
- ◆ 1 phonation valve.
- ◆ 1 cough cap.
- ◆ 1 standard connector.
- ◆ 1 sealing cap.
- ◆ 1 velcro neck band.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 4,0 / 5,0 / 6,0 / 7,0 / 8,0 / 9,0 / 10,0 / 12,0 mm.

CC TRACHEOSTOMY CANNULA WITH CUFF

PRODUCT CODE: 9755

- ◆ 1 PVC cannula, conical shape, curved, polished distal orifice, atraumatic. Pivoting fixation flange. Low-pressure cuff, pilot balloon and retention valve.
- ◆ 2 inner PVC cannulas, interchangeable.
- ◆ 1 insertion aid.
- ◆ 1 cough cap.
- ◆ 1 standard connector.
- ◆ 1 velcro neck band.



STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 5,0 / 5,5 / 6,0 / 6,5 / 7,0 / 7,5 / 8,0 / 8,5 / 9,0 / 9,5 / 10,0 mm.

CC TRACHEOSTOMY CANNULA, FENESTRATED, WITH CUFF (*)

PRODUCT CODE: 9237 (*) - 2 inner cannulas, one with phonation opening and another without phonation opening.

- ◆ 1 PVC cannula, conical shape, curved, longitudinal line opaque to X-Rays, polished distal orifice, atraumatic. Phonation opening. Pivoting fixation flange. Low-pressure cuff, pilot balloon and retention valve.
- ◆ 1 inner PVC cannula without phonation opening.
- ◆ 1 inner PVC cannula with phonation opening.
- ◆ 1 insertion aid.
- ◆ 1 phonation valve.
- ◆ 1 cough cap.
- ◆ 1 standard connector.
- ◆ 1 sealing cap.
- ◆ 1 velcro neck band.



STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 5,0 / 5,5 / 6,0 / 6,5 / 7,0 / 7,5 / 8,0 / 8,5 / 9,0 / 9,5 / 10,0 mm.

CC TRACHEOSTOMY CANNULA, FENESTRATED, WITH CUFF

PRODUCT CODE: 9238



- ◆ 1 PVC cannula, conical shape, curved, longitudinal line opaque to X-Rays, polished distal orifice, atraumatic. Phonation opening. Pivoting fixation flange. Low-pressure cuff, pilot balloon and retention valve.
- ◆ 1 inner PVC cannula.
- ◆ 1 insertion aid.
- ◆ 1 standard connector.
- ◆ 1 cough cap.
- ◆ 1 sealing cap.
- ◆ 1 velcro neck band.

STERILE, SINGLE PACK, SINGLE USE.


AVAILABLE SIZES: 5,0 / 5,5 / 6,0 / 6,5 / 7,0 / 7,5 / 8,0 / 8,5 / 9,0 / 9,5 / 10,0 mm.



CC TRACHEOSTOMY CANNULA WITH CUFF AND SUBGLOTTIC ASPIRATION

PRODUCT CODE: 9259

- ◆ 1 PVC cannula, curved, flexible, transparent, conical, longitudinal line opaque to X-Rays, polished distal orifice, atraumatic. Pivoting fixation flange. Low-pressure cuff, pilot balloon and retention valve. Suction channel with orifice next to the upper edge of the cuff.
- ◆ 2 inner PVC cannulas, interchangeable,
- ◆ 1 insertion aid.
- ◆ 1 cough cap.
- ◆ 1 standard connector.
- ◆ 1 velcro neck band.



STERILE, SINGLE PACK, SINGLE USE.

AVAILABLE SIZES: 5,0 / 5,5 / 6,0 / 6,5 / 7,0 / 7,5 / 8,0 / 8,5 / 9,0 / 9,5 / 10,0 mm.

ADJUSTABLE TRACHEOSTOMY CANNULA WITH CUFF

PRODUCT CODE: 9736



- ◆ 1 PVC cannula, right angle curved, flexible, transparent, longitudinal line opaque to X-rays, polished distal orifice, atraumatic. Standard connector at proximal end. Fixation flange with lock ring, repositionable and adjustable. Low-pressure cuff, pilot balloon and retention valve.
- ◆ 1 suction catheter.
- ◆ 1 velcro neck band.

STERILE, SINGLE PACK, SINGLE USE.

AVAILABLE SIZES: 6,0 / 7,0 / 8,0 / 9,0 / 10,0 / 11,0 mm.

ADJUSTABLE TRACHEOSTOMY CANNULA WITH SUBGLOTTIC ASPIRATION

PRODUCT CODE: 9258

- ◆ 1 PVC cannula, right angle curved, flexible, transparent, longitudinal line opaque to X-rays polished distal orifice, atraumatic. Standard connector at proximal end. Fixation flange with lock ring, repositionable and adjustable. Low-pressure cuff, pilot balloon and retention valve. Suction channel incorporated, with orifice next to the upper edge of the cuff.
- ◆ 1 suction catheter.
- ◆ 1 velcro neck band.

STERILE, SINGLE PACK, SINGLE USE.

AVAILABLE SIZES: 6,0 / 7,0 / 8,0 / 9,0 / 10,0 mm.



TRACHEOSTOMY CANNULA WITH SUBGLOTTIC ASPIRATION

PRODUCT CODE: 9245

- ◆ 1 PVC cannula, curved, flexible, transparent, longitudinal line opaque to X-rays, polished distal orifice, atraumatic. Standard connector at proximal end. Low-pressure cuff, pilot balloon and retention valve. Suction channel, with orifice next to the upper edge of the cuff.
- ◆ 1 insertion aid.
- ◆ 1 velcro neck band.

STERILE, SINGLE PACK, SINGLE USE.

AVAILABLE SIZES: 6,0 / 6,5 / 7,0 / 7,5 / 8,0 / 8,5 / 9,0 / 9,5 / 10,0 mm.



ARMFLEX TRACHEOSTOMY CANNULA

PRODUCT CODE: 9246

- ◆ 1 PVC cannula, reinforced with flexible stainless steel spiral curved, polished distal orifice, atraumatic. Pivoting fixation flange.
- ◆ 2 inner PVC cannulas, interchangeable.
- ◆ 1 insertion aid.
- ◆ 1 phonation valve.
- ◆ 1 cough cap.
- ◆ 1 standard connector.
- ◆ 1 velcro neck band.


STERILE, SINGLE PACK, SINGLE USE.

AVAILABLE SIZES: 7,0 / 8,0 / 9,0 / 10,0 / 11,0 mm.



ARMFLEX FENESTRATED TRACHEOSTOMY CANNULA WITHOUT CUFF

PRODUCT CODE: 9241



- ◆ 1 PVC cannula, reinforced with flexible stainless steel spiral curved, polished distal orifice, atraumatic. Phonation opening. Pivoting fixation flange.
- ◆ 2 inner PVC cannulas, interchangeable, (one with phonation opening and another without phonation opening).
- ◆ 1 insertion aid.
- ◆ 1 cough cap.
- ◆ 1 phonation valve.
- ◆ 1 standard connector.
- ◆ 1 velcro neck band.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 7,0 / 8,0 / 9,0 / 10,0 / 11,0 mm.

ARMFLEX TRACHEOSTOMY CANNULA WITH CUFF

PRODUCT CODE: 9247



- ◆ 1 PVC cannula, reinforced with flexible stainless steel spiral curved, polished distal orifice, atraumatic. Pivoting fixation flange. Low-pressure cuff, pilot balloon and retention valve.
- ◆ 2 inner PVC cannulas, interchangeable,
- ◆ 1 insertion aid.
- ◆ 1 cough cap.
- ◆ 1 standard connector.
- ◆ 1 velcro neck band.


STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 7,0 / 8,0 / 9,0 / 10,0 / 11,0 mm.



ARMFLEX FENESTRATED TRACHEOSTOMY CANNULA WITH CUFF

PRODUCT CODE: 9242

- ◆ PVC cannula, reinforced with flexible stainless steel spiral curved, polished distal orifice, atraumatic. Phonation opening. Pivoting fixation flange. Low-pressure cuff, pilot balloon and retention valve.
- ◆ 2 inner PVC cannulas, interchangeable, (one with phonation opening and another without phonation opening).
- ◆ 1 insertion aid.
- ◆ 1 cough cap.
- ◆ 1 phonation valve.
- ◆ 1 standard connector.
- ◆ 1 velcro neck band.



STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 7,0 / 8,0 / 9,0 / 10,0 / 11,0 mm.

ADJUSTABLE ARMFLEX TRACHEOSTOMY CANNULA WITH CUFF

PRODUCT CODE: 9729



- ◆ 1 PVC cannula reinforced with flexible stainless steel spiral, polished distal orifice, atraumatic. Standard connector at proximal end. Fixation flange with lock ring, repositionable and adjustable. Low-pressure cuff, pilot balloon and retention valve.
- ◆ 1 suction catheter.
- ◆ 1 velcro neck band.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 7,0 / 8,0 / 9,0 / 10,0 / 11,0 mm.



ARMFLEX LARYNGECTOMY TUBE WITH CUFF.

PRODUCT CODE: 9220

- ◆ 1 PVC tube reinforced with flexible stainless steel spiral, polished distal orifice, atraumatic. Low-pressure cuff, pilot balloon and retention valve. Standard connector.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 4,0 / 5,0 / 6,0 / 7,0 / 8,0 / 9,0 / 10,0 mm.



CRICOTOMY SET

PRODUCT CODE: 5980




- ◆ 1 plastic cannula, short and straight, fixation flange, standard connector at proximal end.
- ◆ 1 trocar with steel blade, cutting.
- ◆ 1 tube for connection to ventilation equipment.
- ◆ 1 suction catheter.
- ◆ 1 velcro neck band.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZE: 4,5 mm.



DOUBLE-LUMEN BRONCHIAL TRACHEOSTOMY CANNULA (LEFT - RIGHT)

PRODUCT CODE: 9263 (left L= 9,5 cm) - 9264 (left L= 7,5 cm) - 9265 (left L= 8,5 cm) - 9266 (right L= 8,5 cm) - 9267 (right L= 7,5 cm) - 9268 (right L= 9,5 cm)



- ◆ 1 PVC tube, divided by longitudinal partition in two (2) channels (tracheal and bronchial), polished distal orifices, atraumatic. Semi-seated connectors. Low-pressure cuffs, delimited by radiopaque rings (blue-bronchial / colourless-tracheal), pilot balloon and retention valve. Neck fixation flange.
- ◆ 1 stylet for intubation.
- ◆ 2 angled connection pieces with double sealing cap and standard connector, rotating connectors.
- ◆ 1 plastic "Y" connector.
- ◆ 1 plastic connector extraction key.
- ◆ 1 velcro neck band.
- ◆ Suction catheters with suction control.


STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZE: Ch 39
 Intratracheal limb length: 7,5 / 8,5 / 9,5 cm.

ARMFLEX ENDOTRACHEAL CURVED TUBE WITHOUT CUFF

PRODUCT CODE: 9700

- ◆ 1 PVC curved tube reinforced with flexible stainless steel spiral graduated in cm, polished distal orifice, atraumatic. Murphy eye. Standard connector.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 3,0 / 3,5 / 4,0 / 4,5 / 5,0 / 5,5 / 6,0 / 6,5 / 7,0 / 7,5 / 8,0 / 8,5 / 9,0 / 9,5 / 10,0 mm.



ARMFLEX ENDOTRACHEAL STRAIGHT TUBE WITHOUT CUFF

PRODUCT CODE: 9713

- ◆ 1 PVC straight tube reinforced with flexible stainless steel spiral graduated in cm, polished distal orifice, atraumatic. Standard connector.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 3,0 / 3,5 / 4,0 / 4,5 / 5,0 / 5,5 / 6,0 / 6,5 / 7,0 / 7,5 / 8,0 / 8,5 / 9,0 / 9,5 / 10,0 mm



ARMFLEX ENDOTRACHEAL CURVED TUBE WITH CUFF

PRODUCT CODE: 9703



- ◆ 1 PVC curved tube reinforced with flexible stainless steel spiral graduated in cm, polished distal orifice, atraumatic. Murphy eye. Low-pressure cuff, pilot balloon and retention valve. Standard connector.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 3,0 / 3,5 / 4,0 / 4,5 / 5,0 / 5,5 / 6,0 / 6,5 / 7,0 / 7,5 / 8,0 / 8,5 / 9,0 / 9,5 / 10,0 mm.



ARMFLEX ENDOTRACHEAL STRAIGHT TUBE WITH CUFF

PRODUCT CODE: 9714

- ◆ 1 PVC straight tube reinforced with flexible stainless steel spiral graduated in cm, polished distal orifice, atraumatic. Low-pressure cuff, pilot balloon and retention valve. Standard connector.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 3,0 / 3,5 / 4,0 / 4,5 / 5,0 / 5,5 / 6,0 / 6,5 / 7,0 / 7,5 / 8,0 / 8,5 / 9,0 / 9,5 / 10,0 mm.



ENDOTRACHEAL TUBE WITH / WITHOUT CUFF AND INSTILLATION CHANNEL

PRODUCT CODE: 9217 / 9218



- ◆ 1 PVC tube, flexible, transparent, graduated in cm, longitudinal line opaque to X-rays polished distal orifice, atraumatic. Murphy eye. Low-pressure cuff, pilot balloon and retention valve. Instillation channel in the tube wall with proximal Luer connector and draining orifice distal to cuff. Semi-seated standard connector.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES:
WITH CUFF: 6,5 / 7,0 / 7,5 / 8,0 / 8,5 / 9,0 / 9,5 / 10,0 mm.
WITHOUT CUFF: 2,5 / 3,0 / 3,5 / 4,0 / 4,5 / 5,0 / 5,5 / 6,0 mm.



ENDOTRACHEAL TUBE WITH SUBGLOTTIC ASPIRATION

PRODUCT CODE: 9219




- ◆ 1 PVC tube, flexible, transparent, graduated in cm, longitudinal line opaque to X-rays polished distal orifice, atraumatic. Low-pressure cuff, pilot balloon and retention valve. Murphy eye. Suction channel incorporated, with orifice next to the upper edge of the cuff. Semi-seated standard connector.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 6,0 / 6,5 / 7,0 / 7,5 / 8,0 / 8,5 / 9,0 / 9,5 mm.



ANGLED ENDOTRACHEAL NASAL TUBE WITHOUT CUFF

PRODUCT CODE: 9207




- ◆ 1 PVC tube, anatomic, transparent, graduated in cm, longitudinal line opaque to X-rays polished distal orifice, atraumatic. Murphy eye. Semi-seated standard connector.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 3,0 / 3,5 / 4,0 / 4,5 / 5,0 / 5,5 / 6,0 / 6,5 mm.



ANGLED ENDOTRACHEAL NASAL TUBE WITH CUFF

PRODUCT CODE: 9208




- ◆ 1 PVC tube, anatomic, transparent, graduated in cm, longitudinal line opaque to X-rays polished distal orifice, atraumatic. Murphy eye. Low-pressure cuff, pilot balloon and retention valve. Semi-seated standard connector.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 3,0 / 3,5 / 4,0 / 4,5 / 5,0 / 5,5 / 6,0 / 6,5 / 7,0 / 7,5 / 8,0 / 8,5 mm.



ANGLED ENDOTRACHEAL ORAL TUBE WITHOUT CUFF

PRODUCT CODE: 9205



- ◆ 1 PVC tube, anatomic, transparent, graduated in cm, longitudinal line opaque to X-rays polished distal orifice, atraumatic. Murphy eye.

Semi-seated standard connector.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 3,0 / 3,5 / 4,0 / 4,5 / 5,0 / 5,5 / 6,0 / 6,5 / 7,0 / 7,5 / 8,0 / 8,5 mm.



ANGLED ENDOTRACHEAL ORAL TUBE WITH CUFF

PRODUCT CODE: 9206


- ◆ 1 PVC tube, anatomic, transparent, graduated in cm, longitudinal line opaque to X-rays polished distal orifice, atraumatic. Murphy eye. Low-pressure cuff, pilot balloon and retention valve. Semi-seated standard connector.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 3,0 / 3,5 / 4,0 / 4,5 / 5,0 / 5,5 / 6,0 / 6,5 / 7,0 / 7,5 / 8,0 / 8,5 / 9,0 mm.



STYLETS FOR INTUBATION

PRODUCT CODE: 9710



- ◆ 1 PVC tube with stainless steel core.

NON STERILE, TUBULAR PROTECTIVE CONTAINER, SINGLE USE.


**AVAILABLE SIZES: 1,5 mm x 350 mm
 2,0 mm x 400 mm
 2,6 mm x 400 mm
 3,3 mm x 450 mm
 4,3 mm x 450 mm
 5,6 mm x 500 mm**



DOUBLE-LUMEN BRONCHIAL TUBE (RIGHT)

PRODUCT CODE: 9791 (Ch 33 to 41) - 9783 (Ch 26 and 28)


- ◆ 1 PVC tube, divided by longitudinal partition in two (2) channels (tracheal and bronchial), polished distal orifices, atraumatic. Semi-seated connectors. Low-pressure cuffs, delimited by radiopaque rings (blue-bronchial / colourless-tracheal), pilot balloon and retention valve.
- ◆ 1 stylet for intubation.
- ◆ 2 angled connection pieces with double sealing cap and standard connector, rotating connectors.
- ◆ 1 plastic "Y" connector.
- ◆ 1 plastic connector extraction key.
- ◆ Suction catheter with suction control.



STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: Ch 26 / 28 / 33 / 35 / 37 / 39 / 41.

DOUBLE-LUMEN BRONCHIAL TUBE (LEFT)

PRODUCT CODE: 9792 (Ch 33 to 41) - 9784 (Ch 26 and 28).




- ◆ 1 PVC tube, divided by longitudinal partition in two (2) channels (tracheal and bronchial), polished distal orifices, atraumatic. Semi-seated connectors. Low-pressure cuffs, delimited by radiopaque rings. (blue-bronchial / colourless-tracheal) pilot balloon and retention valve.
- ◆ 1 stylet for intubation.
- ◆ 2 angled connection pieces with double sealing cap and standard connector, rotating connectors.
- ◆ 1 plastic "Y" connector.
- ◆ 1 plastic connector extraction key.
- ◆ Suction catheter with suction control.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: Ch 26 / 28 / 33 / 35 / 37 / 39 / 41.

ESCHMANN INTUBATION STYLETS

PRODUCT CODE: 9230




- ◆ 1 PVC stylet graduated in cm. with curved end.

STERILE, SINGLE PACK.
AVAILABLE SIZES: 2,6 mm x 700 mm
5,2 mm x 700 mm



DOUBLE-LUMEN BRONCHIAL TUBE WITH CARINA HOOK (RIGHT)

PRODUCT CODE: 9772



- ◆ 1 PVC tube, divided by longitudinal partition in two (2) channels (tracheal and bronchial), polished distal orifices, atraumatic. Semi-seated connectors. Low-pressure cuffs, delimited by radiopaque rings (blue-bronchial / colourless-tracheal), pilot balloon and retention valve. Carina hook.
- ◆ 1 stylet for intubation.
- ◆ 2 angled connection pieces with double sealing cap and standard connector, rotating connectors.
- ◆ 1 plastic "Y" connector.
- ◆ 1 plastic connector extraction key.
- ◆ Suction catheter with suction control.


STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: Ch 33 / 35 / 37 / 39 / 41.



DOUBLE-LUMEN BRONCHIAL TUBE WITH CARINA HOOK (LEFT)

PRODUCT CODE: 9773

- ◆ 1 PVC tube, divided by longitudinal partition in two (2) channels (tracheal and bronchial), polished distal orifices, atraumatic. Semi-seated connectors. Low-pressure cuffs, delimited by radiopaque rings (blue-bronchial / colourless-tracheal), pilot balloon and retention valve. Carina hook.
- ◆ 1 stylet for intubation.
- ◆ 2 angled connection pieces with double sealing cap and standard connector, rotating connectors.
- ◆ 1 plastic "Y" connector.
- ◆ 1 plastic connector extraction key.
- ◆ Suction catheter with suction control.




STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: Ch 33 / 35 / 37 / 39 / 41.



URETERAL CATHETER

PRODUCT CODE: 9460 (chevassu) - 9461 (open, curved tip) - 9462 (cylindrical tip, one eye) - 9464 (bevelled tip, one eye) - 9465 (open tip, one eye) - 9466 (open tip, two eyes) - 9467 (open curved tip, one eye) - 9468 (bevelled tip, two eyes) - 9470 (cylindrical tip, two eyes) - 9471 (cylindrical, curved tip, 2 eyes) - 9472 (open tip) - 9473 (olive tip, two eyes) - 9474 (open, curved tip, two eyes) - 9475 (curved, cylindrical tip, one eye)



- ◆ 1 PUR catheter, with wide variety of tips, radiopaque, graduated in cm (for diagnostic and therapeutic procedures).
- ◆ 1 guide wire made of stainless steel.
- ◆ 1 Luer connector.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: Ch 3,0 / 4,0 / 5,0 / 6,0 / 7,0 / 8,0.
Exception: code 9460 size Ch 3,0 not available



DOUBLE "J" URETERAL CATHETER - SET


PRODUCT CODE: 9400 (12 cm)* - 9406 (16 cm)* - 9409 (18 cm)* - 9412 (20 cm)* - 9413 (14 cm)* - 9415 (22 cm)* - 9418 (24 cm)* - 9421 (26 cm)* - 9424 (28 cm)* - 9427 (30 cm)* / 9411 (30 cm)** - 9430 (12 cm)** - 9433 (14 cm)** - 9436 (16 cm)** - 9439 (18 cm)** - 9442 (20 cm)** - 9445 (22 cm)** - 9448 (24 cm)** - 9451 (26 cm)** - 9454 (28 cm)**

- ◆ 1 PUR catheter, radiopaque, graduated in cm, coiled cylindrical tips, black longitudinal positioning line, black circumferential mark at the proximal end, open or closed distal catheter end, drainage orifices in its entirety.
- ◆ 1 guide wire made of stainless steel, PTFE coating.
- ◆ 1 introducer sheath.
- ◆ 2 fixation clamps (for closed tip).

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: Ch 3,0 / 4,8 / 5,0 / 6,0 / 7,0 / 8,0 / 9,0.
Length: 12 / 14 / 16 / 18 / 20 / 22 / 24 / 26 / 28 / 30 cm.



*Closed tip **Open tip



SIMPLE "J" URETERAL CATHETER - SET


PRODUCT CODE: 9457

- ◆ 1 PUR catheter, radiopaque, graduated in cm, coiled cylindrical tip at the renal end, black longitudinal positioning line, drainage orifices at the renal end.
- ◆ 1 guide wire made of stainless steel, PTFE coating.
- ◆ 1 fixation clamp.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: Ch 4,8 / 6,0 / 7,0 / 8,0.
Length: 70 cm.



STONE EXTRACTOR



Code	Wires	Length	Code	Wires	Length
5902	4	75	5943	4	90
5903	5	75	5944	5	90
5904	6	75	5945	6	90

- ◆ 1 Nitinol basket, 4, 5 and 6 wires, flexible helical shape, with atraumatic cylindrical tip.
- ◆ 1 sheath contain shaft and basket.
- ◆ 1 piece right angles, with Luer connector for injection of rinsing solutions or radiopaque contrast liquid.
- ◆ 1 Luer-lock cap.
- ◆ 1 ergonomic handle, slidable, manageable with one hand.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: Ch 3,0 / 4,0 / 5,0.
Length: 75 / 90 cm.



PHthalate
FREE



LATEX
FREE

Code	Wires	Length
5907	4	90

- ◆ 1 Nitinol basket, 4 wires, flexible helical shape, with atraumatic tip.
- ◆ 1 sheath containing shaft and basket.
- ◆ 1 ergonomic handle, slidable, manageable with one hand.


STERIL, SINGLE PACK, SINGLE USE
AVAILABLE SIZES: Ch 2,2 / 2,5
Length: 90 cm

NEPHROSTOMY SET WITH PIGTAIL CATHETER

PRODUCT CODE: 9692

9693 - (Presentation without tray)

- ◆ 1 PUR radiopaque catheter, smooth, flexible, renal end rolled up with drainage orifices, cylindrical coaxial tip with central opening, Support piece for introduction. Detachable Luer-lock connector.
- ◆ 1 percutaneous puncture set
- ◆ 1 guide wire made of stainless steel, PTFE coating (safety "J" tip and piece of rectification).
- ◆ 1 set of dilators (progressive diameters).
- ◆ 1 connection tube with: Luer-lock to the catheter, funnel connection to urine bag.
- ◆ 1 two-way stop-cock.



LATEX
FREE




PHthalate
FREE

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: Ch 6,0 / 8,0 / 10,0 / 12,0.

NEPHROSTOMY SET WITH MALECOT CATHETER

PRODUCT CODE: 9689

9694 - (Presentation without tray)




- ◆ 1 PUR catheter with radiopaque line, smooth, flexible, Malecot type renal end, cylindrical coaxial tip with central opening, Funnel in proximal end. Plastic stylet.
- ◆ 1 percutaneous puncture set.
- ◆ 1 guide wire made of stainless steel, PTFE coating (safety "J" tip and piece of rectification).
- ◆ 1 dilation set (progressive diameters).
- ◆ 1 connection tube with closing clamp progressive tapered end toward the catheter funnel connection to urine bag.

STERILE, SINGLE PACK, SINGLE USE.

AVAILABLE SIZES: Ch 10,0 / 12,0 / 14,0 / 16,0.

CATHETER FOR URETEROCUTANEOSTOMY

PRODUCT CODE: 9401 (26 cm) - 9402 (32 cm)



- ◆ 1 PUR catheter, flexible, multi-fenestrated, with antiperistaltic anchorage fins, and disc for skin fixation.

STERILE, SINGLE PACK, SINGLE USE.

AVAILABLE SIZES: 3,5 / 4,2 mm.

Length: 26 / 32 cm.



CYSTOMETRY CATHETER


PRODUCT CODE: 9504 (adult) - 9503 (pediatric)

- ◆ 1 PVC double-lumen catheter.
Filling lumen with two orifices in the distal end with blue Luer connector.
Lumen for registration with one orifice at distal end with crystal coloured Luer connector.
Introduction marks.

STERILE, SINGLE PACK, SINGLE USE.


AVAILABLE SIZES: Pediatric Ch 6,0 / 7,0 / Length: 32 cm.

Adult Ch 6,0 / 7,0 / 9,0 / Length: 52 cm.



RECTAL TUBE FOR CYSTOMETRY

PRODUCT CODE: 9966



- ◆ 1 PVC tube, crystal color, flexible, Luer connector on the proximal end, Luer-lock cap. Closure clamp. Latex balloon at the distal end.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZE: Ch 12,0 / Length: 38 cm.



FOLEYPLAST, SOFT FOLEYPLAST AND EXTRA SOFT FOLEYPLAST CATHETER

PRODUCT CODE:

9558 - Foleyplast Dufour 30 mL

9533 - Soft Foleyplast Dufour 30 mL

9587 - Extra Soft Foleyplast Dufour 30 mL

9536 - Foleyplast Dufour 75 mL

9535 - Soft Foleyplast Dufour 75 mL

9583 - Extra Soft Foleyplast Dufour 75 mL

9534 - Foleyplast Couvelaire 30 mL

9532 - Soft Foleyplast Couvelaire 30 mL

9586 - Extra Soft Foleyplast Couvelaire 30 mL

9554 - Foleyplast Couvelaire 75 mL

9531 - Soft Foleyplast Couvelaire 75 mL

9584 - Extra Soft Foleyplast Couvelaire 75 mL

9580 - Foleyplast Mercier 30 mL

9581 - Soft Foleyplast Mercier 30 mL

9575 - Foleyplast Mercier 75 mL

9539 - Soft Foleyplast Mercier 75 mL

9585 - Extra Soft Foleyplast Mercier 75 mL




- ◆ 1 PVC tube, 3 way, crystal colour, flexible (for postoperative drainage), range of tips, sizes and location of distal holes (Couvelaire, Dufour, Mercier). Valve for Luer and Luer-lock connection. Latex balloon (30 mL / 75 mL). GREEN longitudinal line = Foleyplast
RED longitudinal line = Soft Foleyplast, Extra Soft Foleyplast

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: Ch 16 / 18 / 20 / 22 / 24.




Dufour



Couvelaire



Mercier



UMBILICAL CATHETER

PRODUCT CODE: 0631



- ◆ 1 PUR catheter, graduated in cm, continuous longitudinal radiopaque line. Polished tip, atraumatic. Luer connector at the proximal end.
- ◆ 1 polyethylene protective packaging.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: Ch 2,5 / 3,0 / 3,5 / 4,0 / 5,0 / 6,0 / 7,0 / 8,0.
Length: 40 cm.




EMBOLECTOMY CATHETER

PRODUCT CODE: 9812

- ◆ 1 PUR catheter, radiopaque, with depth marks (colour coded for size). Latex balloon, attached and tied at both ends. Luer connector at the proximal end.
- ◆ 1 stainless steel filament attached to a Luer-lock cap.
- ◆ 1 syringe three parts for insufflation.
- ◆ 1 polyethylene protective packaging.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: Ch 2,0 / 3,0 / 4,0 / 5,0 / 6,0 / 7,0 / 8,0.



TUBE FOR TREATMENT OF ESOPHAGEAL VARICES

PRODUCT CODE: 9825 - 9827 (pediatric)



- ◆ 1 PVC tube, 3-way, rounded distal end, with multiple openings to irrigation and suction.
Gastric latex balloon in connection with insufflation valve with "Stomach" inscription.
Esophageal latex balloon in connection with insufflation valve with "Esophagus" inscription.
Radiopaque rings at the proximal and distal end of both balloons.
Funnel in proximal end.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: Ch 14,0 / 16,0 / 18,0 / 21,0.
Ch 14,0 (pediatric).



RECTAL TUBE WITH PVC BALLOON

PRODUCT CODE: 9961

- ◆ 1 PVC tube, rounded distal end,
side orifice for irrigation and suction.
PVC balloon.
Funnel in proximal end.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 8,0 / 10,0 mm.



RECTAL TUBE WITH LATEX BALLOON

PRODUCT CODE: 9962

- ◆ 1 PVC tube, rounded distal end, side orifice
for irrigation and suction.
Latex balloon.
Funnel in proximal end.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: 8,0 / 10,0 mm.





EN ISO 13485
CE 0124 - Annex II, Section 3 of the Directive 93/42/EEC
Certified since 1998 by DEKRA Certification GmbH - Germany

Servando Gómez 2450
ZIP CODE: 12100 (Montevideo - Uruguay)

Phone: +598 26005211 / +598 26006073 / +598 26005835
FAX: +598 26013958
Email: electroplast@epsa.com.uy

Please visit our website at:
epsa.uy