

**EPSA<sup>®</sup>** *safety first*


CRICOTOMY SET

# CRICOTOMY SET

**CODE: 5980**

- 1 Plastic cannula, short and straight, fixation flange.
- Standard connector at proximal end.
- 1 trocar with steel blade, cutting
- 1 tube for connection to ventilation equipment
- 1 suction catheter
- 1 velcro neck band

**STERILE, SINGLE PACK, SINGLE USE.**



**Available Sizes**

**4,5 mm.**


## Recommendation for

### use

Recommendation of a professional on the use of the technology.

The Cricotomy Set is indicated for rapid access to the airway when other maneuvers have failed or cannot be performed:

- . Sudden complete airway obstruction.
- . Foreign body airway obstruction refractory to attempts to elimination.
- . Complete airway occlusion: by mass lesion (tumor, hematoma and abscess).
- . Severe upper respiratory edema:
  - Anaphylaxis.
  - Thermal inhalation injuries.
  - Caustic ingestion.
  - Angioedema.



## Recommendation for use

Epiglottitis complicated by severe respiratory compromise and / or respiratory arrest.

- . Inability to perform intubation, due to:
  - Bleeding in the upper respiratory tract
  - Anatomical variants (acquired or congenital).
  - Regurgitation and / or massive aspiration.
  - Severe maxillofacial trauma (fractures of the middle part of the face).

