



**EPSA<sup>®</sup>** *safety first*

**FOLEYPLAST, SOFT FOLEYPLAST  
AND EXTRA SOFT FOLEYPLAST CATHETER**

# FOLEYPLAST, SOFT FOLEYPLAST AND EXTRA SOFT FOLEYPLAST CATHETER

1 PVC tube, 3-way, crystal colour, flexible,  
(for postoperative drainage) range of tips, sizes and  
location of distal holes.

(Dufour, Couvelaire, Mercier)


Valve for Luer and Luer-lock connection.

Latex balloon (30 ml. / 75 ml.)

GREEN longitudinal line = FoleyPlast

RED longitudinal line = SOFT FoleyPlast, Extra Soft  
Foleyplast

**STERILE, SINGLE PACK, SINGLE USE.**



**Available Sizes**

Ch 16,0 / 18,0 / 20 ,0 / 22,0 / 24,0

# FOLEYPLAST, SOFT FOLEYPLAST AND EXTRA SOFT FOLEYPLAST CATHETER

9558	Foleyplast Dufour 30 mL	9532	Soft Foleyplast Couvelaire 30 mL
9533	Soft Foleyplast Dufour 30 mL	9554	Foleyplast Couvelaire 75 mL
9536	Foleyplast Dufour 75 mL	9531	Soft Foleyplast Couvelaire 75 mL
9535	Soft Foleyplast Dufour 75 mL	9580	Foleyplast Mercier 30 mL
9534	Foleyplast Couvelaire 30 mL	9581	Soft Foleyplast Mercier 30 mL
9575	Foleyplast Mercier 75 mL	9539	Soft Foleyplast Mercier 75 mL
9583	Extra Soft Foleyplast Dufour 75mL	9584	Extra Soft Foleyplast Couvelaire 75mL
9585	Extra Soft Foleyplast Mercier 75mL		

## Recommendation for use

Recommendation of a professional on the use of the technology.

The FoleyPlast, Soft FoleyPlast and Extra Soft FoleyPlast Catheter are indicated in the postoperative period of:

- \* Bladder surgeries, such as:
  - partial cystectomy, – resection of tumors.
  - fulguration of tumors,
- \* RTU
- \* Injuries associated with: – gynecological surgeries – other urological procedures.

