

EPSA[®] *safety first*

**NEPHROSTOMY SET
WITH MALECOT CATHETER**

NEPHROSTOMY SET WITH MALECOT CATHETER

CODE: 9689

1 PUR catheter with radiopaque line, smooth, flexible, Malecot tip renal end, cylindrical coaxial tip with central opening. Funnel in proximal end. Plastic stylet.

1 percutaneous puncture set.

1 guide wire made of stainless steel, PTFE coating (safety "J" tip and piece of rectification).

1 dilation set (progressive diameters).

1 connection tube with closing clamp progressive tapered end toward the catheter funnel connection to urine bag.



Available Sizes

Ch 10,0 / 12,0 / 14,0 / 16,0

STERILE, SINGLE PACK, SINGLE USE.

Recommendation for use

Recommendation of a professional on the use of the technology.

The Nephrostomy Set with Malecot Catheter is indicated when a nephrostomy is performed, which consists of communicating the urinary excretory tract with the outside, at the level of the renal cavities, when there is an alteration or interruption of the trans-ureteral urinary flow. The application of the nephrostomy catheter allows:

- decompress an obstructed system and maintain or improve kidney function impaired by

ureteral obstruction,

- obtain access to the renal pelvis for radiological procedures.

