



EPSA[®] *safety first*


**NEPHROSTOMY SET
WITH PIGTAIL CATHETER**

NEPHROSTOMY SET WITH PIGTAIL

CODE: 9692

- 1 Nitinol basket, 4 wires, flexible helical shape, with atraumatic tip.
- 1 sheath containing shaft and basket.
- 1 ergonomic handle, slidable, manageable with one hand.

STERILE, SINGLE PACK, SINGLE USE.



Available Sizes
Ch 6,0 / 8,0 / 10,0 / 12,0

Recommendation for use

Recommendation of a professional on the use of the technology.

The Nephrostomy Set with Pigtail Catheter is indicated when a nephrostomy is performed, which consists of communicating the urinary excretory tract with the outside, at the level of the renal cavities, when there is an alteration or interruption of the trans-ureteral urinary flow. The application of the nephrostomy catheter allows:

- decompress an obstructed system and maintain or improve kidney function impaired by

ureteral obstruction,

- obtain access to the renal pelvis for radiological procedures.

