

**EPSA<sup>®</sup>** *safety first*

**PEDIATRIC AND NEONATAL  
TRACHEOSTOMY CANNULA  
WITHOUT CUFF**

# PEDIATRIC AND NEONATAL TRACHEOSTOMY CANNULA WITHOUT CUFF

CODE: 9764 / 9758

## PEDIATRIC

1 PVC cannula, curved, flexible, transparent, longitudinal line opaque to X-rays, polished tip, atraumatic. Flexible fixation flange with standard connector

1 cloth tape for neck fixation

STERILE, SINGLE PACK, SINGLE USE

## NEONATAL

1 PVC cannula, curved, flexible, transparent, longitudinal line opaque to X-rays, polished tip, atraumatic. Flexible fixation flange with standard connector.

1 cloth tape for neck fixation

**STERILE, SINGLE PACK, SINGLE USE.**



## Pediatric Sizes

3,0 / 3,5 / 4,0 / 4,5 mm

## Neonatal Sizes

2,5 / 3,0 / 3,5 / 4,0 / 4,5 mm

## Recommendation for use

Recommendation of a professional on the use of the technology.

The Pediatric and Neonatal Tracheostomy Cannula without cuff is indicated in tracheostomized patients to allow the maintenance of respiratory function, with adequate management of the airway, in relation to hydration, aspiration of secretions and endoscopic maneuvers for diagnosis or treatment.

It is supplied in a special acrylic container for protection during storage.

