



EPSA[®] *safety first*

**RECTAL TUBE
FOR CYSTOMETRY**

RECTAL TUBE FOR CYSTOMETRY

CODE: 9966

1 PVC tube, cristal color, flexible. Luer connector on the proximal end. Luer-lock cap. Closure clamp. Latex balloon at the distal end.

STERILE, SINGLE PACK, SINGLE USE.



Available Sizes
Ch 12,00 / 38,00 cm

Recommendation for use

Recommendation of a professional on the use of the technology.

The Rectal Tube for Cystometry is indicated for measuring rectal pressure during urodynamic studies.

It allows: to record intra-abdominal pressure variations during the development of the study.

