

EPSA[®] *safety first*


**RECTAL TUBE WITH PVC
BALLOON**

RECTAL TUBE WITH PVC BALLOON

CÓDIGO: 9961

1 PVC tube, rounded distal end, side orifice for irrigation and suction. PVC balloon. Funnel in proximal end.

STERILE, SINGLE PACK, SINGLE USE



AVAILABLE SIZES

8,0 / 10,0 mm

Recomendaciones de uso

The Rectal Tube with PVC Balloon is indicated for diagnostic and therapeutic purposes in order to:

- Administration of enemas: cleaning and preparation for imaging studies; medicated and antiseptic.
- Administration of contrast media to carry out radiological studies.
- Treatment of intestinal obstruction with abdominal distention.
- Reduction of body temperature.

