

EPSA[®] *safety first*

**SIMPLE “J” URETERAL
CATHETER – SET**

SIMPLE “J” URETERAL CATHETER – SET

CODES: 9457

1 PUR catheter, radiopaque, graduated in cm., coiled cylindrical tip at the renal end, black longitudinal positioning line, drainage orifices at the renal end
1 guide wire made of stainless steel, PTFE coating
1 fixation clamp

STERILE, SINGLE PACK, SINGLE USE.



Available Sizes
Ch 4,8 / 6,0 / 7,0 / 8,0.
Length: 70 cm

Recommendation for use

Recommendation of a professional on the use of the technology.


The Simple “J” Ureteral Catheter Set is indicated:

– To guarantee the drainage of the renal collecting system and maintain a guardian of the

ureter when the obstruction or area of the scar is in that sector of the urinary system.

It can be used after performing:

- endopyelotomy
- endoscopic percutaneous ureterotomy
- ureterolithotomy of impacted stones in the upper part of the ureter.
- In reconstructive surgery after radical cystectomy.



Various types of reconstructive surgery can be done, depending on the medical condition

and personal preferences:

- Incontinent referral
- Continent bypass
- Neo-bladder