

**EPSA<sup>®</sup>** *safety first*


**STONE EXTRACTOR**  
– CH 2,2 – 2,5

# STONE EXTRACTOR – CH 2,2 – 2,5

CODE	WIRES	LENGTH
5907	4	90

- 1 Nitinol basket, 4 wires, flexible helical shape, with atraumatic tip.
- 1 sheath containing shaft and basket.
- 1 ergonomic handle, slidable, manageable with one hand.

**STERILE, SINGLE PACK, SINGLE USE.**



**Available Sizes**

Ch 2,2 / 2,2 / 2,5.

Longitud: 90 cm.

## Recommendation for use

Recommendation of a professional on the use of the technology.

The Stone Extractor is indicated for the endoscopic extraction of stones in canalicular organs, namely:

- urinary tract (renal cavities, ureter, bladder, urethra).
- biliary (hepatic ducts, common bile duct, cystic duct, gallbladder).
- pancreatic excretory duct. – excretory ducts of salivary glands.

