

EPSA[®] *safety first*


STONE EXTRACTOR
– CH 1,8 – 2,5

STONE EXTRACTOR - CH 1,8 - 2,5

| CODE | WIRES | LENGTH |
|------|-------|--------|
| 5907 | 4 | 90 |

- 1 Nitinol basket, 4 wires, flexible helical shape, with atraumatic tip.
- 1 sheath containing shaft and basket.
- 1 ergonomic handle, slidable, manageable with one hand.

STERILE, SINGLE PACK, SINGLE USE.



Available Sizes

Ch 1,8 / 2,2 / 2,5.

Longitud: 90 cm.

Recommendation for use

Recommendation of a professional on the use of the technology.

The Stone Extractor is indicated for the endoscopic extraction of stones in canalicular organs, namely:

- urinary tract (renal cavities, ureter, bladder, urethra).
- biliary (hepatic ducts, common bile duct, cystic duct, gallbladder).
- pancreatic excretory duct. – excretory ducts of salivary glands.

