

EPSA[®] *safety first*


STONE EXTRACTOR
– CH 3,0 – 5,0

STONE EXTRACTOR – CH 3,0 – 5,0

- 1 Nitinol basket 3, 4, 5 and 6 wires, flexible helical shape, with atraumatic cylindrical tip.
- 1 sheath contain shaft and basket
- 1 piece right angles, with Luer connector for injection of rinsing solutions or radiopaque contrast liquid
- 1 Luer-lock cap
- 1 ergonomic handle, slidable, manageable with one hand

STERILE, SINGLE PACK, SINGLE USE.

CODE	WIRES	LENGTH	CODE	WIRES	LENGTH
5902	4	75	5940	4	110
5903	5	75	5942	6	110
5904	6	75	5943	4	90
5939	5	110	5944	5	90
			5945	6	90



Available Sizes

Ch 3,0 / 4,0 / 5,0.

Length: 75 / 90 / 110 cm.

Recommendation for use

The Stone Extractor is indicated for the endoscopic extraction of stones in canalicular organs, namely:

- urinary tract (renal cavities, ureter, bladder, urethra).
- biliary (hepatic ducts, common bile duct, cystic duct, gallbladder).
- pancreatic excretory duct. – excretory ducts of salivary glands.

