




EPSA[®] *safety first*

**RECTAL TUBE WITH
LATEX BALLOON**

RECTAL TUBE WITH LATEX BALLOON

CODE: 9962

1 PVC tube, rounded distal end, side orifice for irrigation and suction. Latex balloon. Funnel in proximal end.



Available Sizes

8,0 / 10,0 mm

STERILE, SINGLE PACK, SINGLE USE.

Recommendation for use

The Rectal Tube with Látex Balloon is indicated for diagnostic and therapeutic purposes in order to:

- Administration of enemas: cleaning and preparation for imaging studies; medicated and antiseptic.
- Administration of contrast media to carry out radiological studies.
- Treatment of intestinal obstruction with abdominal distention.
- Reduction of body temperature.

